



## FE Focus B II

### Instruction

【Product name】 FE Focus B II

【Package specification】 1×20 mL, 1×22 mL, 1×25 mL, 1×60 mL, 1×125 mL.

【Intended use】

This product is used to determine position of focal plane of microscopic imaging system for Urine Sediment Analyzer, FUS Series Urinalysis Hybrid and MUS Series Urinalysis System.

【Test principle】 Flow cytometry imaging

【Main component】 QC blood: 2,500 particles/μL~3,500 particles/μL;

**Note:** fixed values of different batches of formed elements analysis focus are different (see the label for details).

【Storage conditions and shelf life】

1.This product shall be sealed and kept in a dark place at the temperature of 2°C~8°C. For shelf life refer to the label.

2.After being unsealed for the first time, it shall be sealed and kept in a dark place at the temperature of 2°C~8°C, and the open-vial shelf life is 30 days.

【Test method】 Refer to the user manual of the applicable instrument.

【Product performance indices】

1. Accuracy: total particles: 2,500 particles/μL~3,500 particles/μL;

2. Homogeneity: CV≤15.0%

【Cautions】

1. This product is used for in vitro diagnostics only;

2. Shake the product for 1~3min before use, and do not shake it violently. If it is not being mixed well, the tested number of particles may be smaller, which will affect the focus effect.

3. Freeze is prohibited.

4.To avoid or reduce relevant infection risks, please wear disposable gloves and eyes/face protection articles when using the product.

【References】

1.Ferris JA urine microscopy comparison & standard: Lab Med14: 659-662,1983.

2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

【Instruction Approved and Modified Date】 05/2023

【Symbol Explanation】



Consult instructions for use



Temperature limit



Manufacturer



In vitro diagnostic medical device



Batch code



Use-by date



Keep away from sunlight



Authorised Representative



Catalogue Number



Comply with In Vitro Diagnostic Devices Directive (98/79/EC)



This Way Up

### DIRUI INDUSTRIAL CO.,LTD.

95 Yunhe Street, New & High Tech.

Development Zone Changchun, Jilin 130012 P.R.China

Tel: +86 431 85100409

Fax: +86 431 85172581

E-mail: dirui@dirui.com.cn <http://www.dirui.com.cn>



Emergo Europe B.V. Westervoortsedijk 60  
6827 AT Arnhem The Netherlands