MAS[®] UA DipTube

LIQUID ASSAYED URINALYSIS CONTROL

IVD

INTENDED USE

MAS[®] UA DipTube is intended for use in the clinical laboratory as a control for qualitative and semi-quantitative procedures used in routine urinalysis testing. Assay values are provided for specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and equipment. The product is suitable for use as a control material for physiochemical, chemical, and microscopic methods of routine urinalysis. **UA DipTube** may be used in conjunction with commercially available standardized systems for urine microscopic analysis.

PRODUCT DESCRIPTION

UA DipTube is a liquid stable control material prepared from human urine. Analyte levels are adjusted with various pure chemicals and human source materials. **UA DipTube** also contains preservatives and stabilizers.

CAUTION: UA DipTube is prepared from human urine pools and contains other human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Antibody (HCV) and HIV-1 and HIV-2 Antibody. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from the Center for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories," 2009.

STORAGE AND STABILITY

Refrigerated Storage (2-8°C): Store unopened **UA DipTube** at 2-8°C. Unopened tubes are stable until the expiration date on the label. Once opened, tubes of control are stable for 8 weeks or 20 reagent strip immersions, whichever occurs first.

Room Temperature Storage (18-25°C): Open tubes of control are stable 2 weeks or 10 reagent strip immersions, whichever comes first.

Do not freeze.

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

EXPECTED VALUES

Expected values may vary slightly with different lot numbers and reagent strips. Refer to the included table for values obtained on specific systems. Good laboratory practice suggests that each laboratory establish its own parameters.

Values listed were obtained using instrument(s) and instrument manufacturer's reagents and procedures at the time of assay. Subsequent modifications in instrument, reagent, or procedures may invalidate these values.

Laboratories using methods other than those listed should establish their own expected ranges.

INSTRUCTIONS FOR USE

Materials Required But NOT Provided:

- Slides
- Pipettors
- Centrifuge tubes, slides, and pipettors from standardized urinalysis systems may be used.
- Centrifuge capable of generating an RCF of 400 for KOVA method or 450 for manual methods.
- Microscope with a 10x objective and a 40x objective lens
- Urinalysis stain (Optional)
- Chemical reagent strips
- Any additional laboratory equipment or reagents required to perform the physiochemical, chemical and microscopic analyses required.

Allow the **UA DipTube** to come to room temperature. Thoroughly mix the contents of the tube before each use by gently inverting for several minutes (3 minutes minimum). Replace cap immediately and store the opened tube at 2-8°C or at room temperature (18-25°C).

Physical Tests:

1. **Appearance:** Visually inspect the color and turbidity.
2. **Specific Gravity:** Measure the specific gravity using a refractometer or reagent strip.

Chemical Tests:

1. If reading results visually, compare the reagent strip with the manufacturer's color comparison chart. If using an automated strip reader, follow manufacturer's directions for use.
2. **Reagent Strips:** Remove the cap from the tube and immerse the reagent strip into the control. Remove strip promptly and follow the manufacturer's instructions for testing. Replace cap immediately after use.
3. **Reagent Tablets or Pregnancy Cassettes:** Remove the cap from the tube and pipette the required amount of control per reagent manufacturer's instructions. Replace cap immediately after use.

Microscopic Analysis:

If using a standardized methodology for microscopic analysis, follow manufacturer's recommendations. Each laboratory should establish its own values for systems other than those listed.

If using a non-standardized methodology for microscopic analysis, centrifuge the **UA DipTube** for 5 minutes at 450 RCF (relative centrifugal force). Pour off supernatant leaving approximately 0.5 mL. Gently resuspend sediment transferring one drop onto a clean, dry microscope slide and cover with a 22 mm coverslip. Values reported are the range of microscopic elements counted in 10 hpf (high powered fields). Each laboratory should establish its own ranges.

Handle the **UA DipTube** in the same manner as patient specimens. Consistent handling of the control is necessary for uniform results.

This product contains cellular elements which may form a small amount of sediment in the bottle. To ensure consistent recoveries, invert and mix well prior to use.

The tube contains 12mL of control. **A fresh tube or a tube that has had no more than 1 reagent strip immersion should be used for microscopic analysis.** A consistent volume for testing is essential for reproducible microscopic results.

QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

LIMITATIONS OF PROCEDURE

UA DipTube must be mixed well before each use to ensure reproducible results.

A consistent volume must be used for testing to ensure reproducible microscopic results. A fresh tube or a tube that has had no more than 1 reagent strip immersion should be used for microscopic analysis.

Organic spheres used to simulate leukocytes may not exhibit the same staining characteristics of human leukocytes.

See manufacturer's product insert for a listing of limitations for each test method. For professional use only.

TECHNICAL ASSISTANCE

In the USA, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit www.maslablink.com and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

Outside of the USA, if your laboratory subscribes to LabLink xL, visit www.maslablink.com select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
UAT-MP	UA DipTube, Multi-Pack	10 x 12 mL

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For Insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For Insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For Insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
 LVL 2 2017-02-28

Thermo
 SCIENTIFIC

For Insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For Insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
LVL 2 2017-02-28

Thermo
SCIENTIFIC

For Insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
 LVL 2 2017-02-28

Thermo
 SCIENTIFIC

For Insert updates go to:
www.thermoscientific.com/diagnostics

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LVL 1 2017-02-28
 LVL 2 2017-02-28

Thermo
 SCIENTIFIC

For Insert updates go to:
www.thermoscientific.com/diagnostics

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IVD

 (a)			
		LVL 1 ^(a) <R> UAT17021	LVL 2 ^(b) <R> UAT17022
KOVA (Qualitative)	Red Cells	0 hpf	1 - 62 hpf
	White Cells	0 hpf	1 - 57 hpf
(Quantitative)	Crystals	Absent	Absent - Present
	Casts	0 ^(c)	0 ^(c)
Cen-Slide	Red Cells	0 hpf	1 - 65 hpf
	White Cells	0 hpf	1 - 64 hpf
Yellow IRIS	Crystals	Absent	Absent - Present
	Casts	0 ^(c)	0 ^(c)
IRIS IQ200	Red Cells	**	**
	White Cells	**	**
URI-system	Crystals	**	**
	Casts	**	**
Coverslide Method	Red Cells	0 - 5 hpf	6 - >100 hpf
	White Cells	0 - 5 hpf	3 - 75 hpf
Yellow IRIS	Crystals	Absent	Absent - Present
	Casts	0 ^(c)	0 ^(c)
URI-system	Red Cells	0 hpf	1 - 90 hpf
	White Cells	0 hpf	1 - 85 hpf
Coverslide Method	Crystals	Absent	Absent - Present
	Casts	0 ^(c)	0 ^(c)

Acetest®, Advantus, Atlas, Clinitek®, Series, Clinitest®, Ictotest®, Micro-Bumintest®, Status®, Diastix®, Reg. TM: Siemens Healthcare Diagnostic, Tarrytown, NY
Array 360®, ICON®, Image®, Reg. TM: Beckman Coulter Inc., Fullerton, CA
Biovation, Reg. TM: Biovation, Inc., Hercules, CA
Dimension®, Reg. TM: Siemens Healthcare Diagnostic, Glasgow, DE
URI-system®, Sure-View, Reg. TM: Thermo Fisher Scientific, Waltham, MA
Signify®, Contrast, OSOM, Reg. TM: Genzyme Corporation, Cambridge, MA
KOVA®, Reg. TM: Hycor Biomedical Inc., Garden Grove, CA
DiaScreen Series, Reg. TM: Hypoguard, Minneapolis, MN
Iris IQ200, Yellow IRIS®, Reg. TM: International Remote Imaging Systems, Inc., Chatsworth, CA
Chemstrip®, Chemstrip® 101, Criterion, CUA, Hitachi®, Micral, Mini UA, SUA, URISYS 1100, URISYS Systems, Reg. TM: Roche Diagnostics, Indianapolis, IN
QuickVue®, UrinQuick, Reg. TM: Quidel Corporation, San Diego, CA
Cen-Slide®, Reg. TM: StatSpin, Norwood, MA
PocketChem UA®, Reg. TM: ThermoFisher Scientific, Boulder, CO
YD Uriscan®, Reg. TM: YD Diagnostics, Seoul, Korea

* **Atypical colors may be observed. Results based on intensity.**
* Atypische Farben können beobachtet werden. Resultate basiert worden auf Intensität.
* On peut observer des couleurs atypiques. Résultats basés sur l'intensité.
* Si possono osservare colori atipici. Risultati basati sull'intensità.
* Pueden observarse colores anormales. Los resultados se basan en la intensidad.
* Der kan observeres atypiske farver. Resultater baseret på intensitet.
* Er kan sprake zijn van atypische kleuren. De resultaten zijn gebaseerd op de intensiteit.
* Epätyypillisiä värejä voidaan havaita. Tulokset perustuvat voimakkuuteen.
* Ενδέχεται να παρατηρηθούν άτυπα χρώματα. Τα αποτελέσματα βασίζονται στην ένταση.
* Det kan forekomme atypiske farger. Resultatene er basert på intensitet.
* Może wystąpić niespodziewane zabarwienie. Wyniki na podstawie intensywności.
* Podem ser observadas cores atípicas. Resultados com base na intensidade da cor.
* Наблюдаемый цвет может быть необычным. Результаты зависят от интенсивности.
* Atypiska färger kan observeras. Resultat baserade på intensitet.
* Atipik renkler gözlemlenebilir. Sonuçlar yoğunluğa dayanır.
* 可能观察到非典型颜色。结果依据强度而定。

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


If you have any questions concerning this notification, in the USA please contact the Technical Support Department at 800-232-3342 or 510-979-5417. Outside the USA, please contact your local Subsidiary or Distributor.

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USA
☎ 800-232-3342
sales.diagnostics.fmt@thermofisher.com

Canada
☎ 800-282-4075
info.cddcanada@thermofisher.com

Germany
☎ +49 0800-40 40 771
cdx.de.info@thermofisher.com

Spain, Portugal & Italy
☎ +34 93589 8338
cdx.es.info@thermofisher.com

Asia Pacific
☎ +61 1800 333 110
cdd.asia.info@thermofisher.com

China
☎ +86 800-810-5118
cdx.cn.info@thermofisher.com

New Zealand
☎ 0800 933 966
auiinfo@thermofisher.com

Switzerland & Austria
☎ +41 26 663 86 70
cdx.ch.info@thermofisher.com

Australia
☎ +61 1800 333 110
auiinfo@thermofisher.com

France
☎ +33 1 40 86 65 20
cdx.fr.info@thermofisher.com

Nordic
☎ +47 2 325 0433
info.nordic.cdd@thermofisher.com

United Kingdom & Ireland
☎ +44 1442 868 940
cdx.uk.info@thermofisher.com

For countries not listed:
distributor.cdd@thermofisher.com

Microgenics Corporation
46500 Kato Road
Fremont, CA 94538 USA
ISO 13485 Certified Company
www.thermoscientific.com/diagnostics

EC REP

Microgenics GmbH
Spitalhofstrasse 94
D-94032 Passau, Germany
Tel: +49 (0) 851 886 89 0
Fax: +49 (0) 851 886 89 10

UAT-INS-VA
Rev. 11 2014 08

Thermo
SCIENTIFIC