# MAS<sup>®</sup> UA DipTube

## LIQUID ASSAYED URINALYSIS CONTROL

IVD

### INTENDED USE

**MAS<sup>®</sup> UA DipTube** is intended for use in the clinical laboratory as a control for qualitative and semi-quantitative procedures used in routine urinalysis testing. Assay values are provided for specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and equipment. The product is suitable for use as a control material for physicochemical, chemical, and microscopic methods of routine urinalysis. **UA DipTube** may be used in conjunction with commercially available standardized systems for urine microscopic analysis.

### PRODUCT DESCRIPTION

**UA DipTube** is a liquid stable control material prepared from human urine. Analyte levels are adjusted with various pure chemicals and human source materials. UA DipTube also contains preservatives and stabilizers.

**DANGER: UA DipTube** control contains <0.9% bovine serum albumin (BSA), ≤0.1% Sulfbromophthalein sodium salt and ≤0.006% 3:1 Mixture: 5-chloro-2-methyl-4-isothiazolin-3-one + 2-methyl-4-isothiazolin-3-one (PROCLIN).

H317 - May cause allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection/face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

**UA DipTube** is prepared from human urine pools and contains other human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Antibody (HCV) and HIV-1 and HIV-2 Antibody. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from the Center for Disease Control/ National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories," 2009.

### STORAGE AND STABILITY

**Refrigerated Storage (2-8°C):** Store unopened **UA DipTube** at 2-8°C. Unopened tubes are stable until the expiration date on the label. Once opened, tubes of control are stable for 8 weeks or 20 reagent strip immersions, whichever occurs first.

**Room Temperature Storage (18-25°C):** Open tubes of control are stable 2 weeks or 10 reagent strip immersions, whichever comes first.

#### Do not freeze.

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

### EXPECTED VALUES

Expected values may vary slightly with different lot numbers and reagent strips. Refer to the included table for values obtained on specific systems. Good laboratory practice suggests that each laboratory establish its own parameters.

Values listed were obtained using instrument(s) and instrument manufacturer's reagents and procedures at the time of assay. Subsequent modifications in instrument, reagent, or procedures may invalidate these values.

Laboratories using methods other than those listed should establish their own expected ranges.

### INSTRUCTIONS FOR USE

#### Materials Required But NOT Provided:

- Slides
- Pipettors
- Centrifuge tubes, slides, and pipettors from standardized urinalysis systems may be used.
- Centrifuge capable of generating an RCF of 400 for KOVA method or 450 for manual methods.
- Microscope with a 10x objective and a 40x objective lens
- Urinalysis stain (Optional)
- Chemical reagent strips
- Any additional laboratory equipment or reagents required to perform the physicochemical, chemical and microscopic analyses required.

Allow the **UA DipTube** to come to room temperature. Thoroughly mix the contents of the tube before each use by gently inverting for several minutes (3 minutes minimum). Replace cap immediately and store the opened tube at 2-8°C or at room temperature (18-25°C).

#### Physical Tests:

1. **Appearance:** Visually inspect the color and turbidity.
2. **Specific Gravity:** Measure the specific gravity using a refractometer or reagent strip.

#### Chemical Tests:

1. If reading results visually, compare the reagent strip with the manufacturer's color comparison chart. If using an automated strip reader, follow manufacturer's directions for use.
2. **Reagent Strips:** Remove the cap from the tube and immerse the reagent strip into the control. Remove strip promptly and follow the manufacturer's instructions for testing. Replace cap immediately after use.
3. **Reagent Tablets or Pregnancy Cassettes:** Remove the cap from the tube and pipette the required amount of control per reagent manufacturer's instructions. Replace cap immediately after use.

#### Microscopic Analysis:

If using a standardized methodology for microscopic analysis, follow manufacturer's recommendations. Each laboratory should establish its own values for systems other than those listed.

If using a non-standardized methodology for microscopic analysis, centrifuge the UA DipTube for 5 minutes at 450 RCF (relative centrifugal force). Pour off supernatant leaving approximately 0.5 mL. Gently resuspend sediment transferring one drop onto a clean, dry microscope slide and cover with a 22 mm coverslip. Values reported are the range of microscopic elements counted in 10 hpf (high powered fields). Each laboratory should establish its own ranges.

Handle the **UA DipTube** in the same manner as patient specimens. Consistent handling of the control is necessary for uniform results.

This product contains cellular elements which may form a small amount of sediment in the bottle. To ensure consistent recoveries, invert and mix well prior to use.

The tube contains 12mL of control. **A fresh tube or a tube that has had no more than 1 reagent strip immersion should be used for microscopic analysis.** A consistent volume for testing is essential for reproducible microscopic results.

### QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

### LIMITATIONS OF PROCEDURE

**UA DipTube** must be mixed well before each use to ensure reproducible results.

A consistent volume must be used for testing to ensure reproducible microscopic results. A fresh tube or a tube that has had no more than 1 reagent strip immersion should be used for microscopic analysis.

Organic spheres used to simulate leukocytes may not exhibit the same staining characteristics of human leukocytes.

See manufacturer's product insert for a listing of limitations for each test method. For professional use only.

### TECHNICAL ASSISTANCE

In the USA, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit [www.maslablink.com](http://www.maslablink.com) and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

Outside of the USA, if your laboratory subscribes to LabLink xL, visit [www.maslablink.com](http://www.maslablink.com) select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
UAT-MP	UA DipTube, Multi-Pack	10 x 12 mL

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
LVL 2 2019-02-28

**Thermo**  
SCIENTIFIC

For Insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For Insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
LVL 2 2019-02-28

**Thermo**  
SCIENTIFIC

For Insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
LVL 2 2019-02-28

**Thermo**  
SCIENTIFIC

For Insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For Insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
 LVL 2 2019-02-28

**Thermo**  
 SCIENTIFIC

For Insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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LVL 1 2019-02-28  
LVL 2 2019-02-28

**Thermo**  
SCIENTIFIC

For insert updates go to:  
[www.thermoscientific.com/diagnostics](http://www.thermoscientific.com/diagnostics)

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Acetest®, Advantus, Atlas, Clinitek® Series, Clinitest®, Ictotest®, Micro-Bumintest®, Status®, Diastix®, Reg. TM: Siemens Healthcare Diagnostic, Tarrytown, NY  
 Array 360®, ICON®, Image®, Reg. TM: Beckman Coulter Inc., Fullerton, CA  
 Biovation, Reg. TM: Biovation, Inc., Hercules, CA  
 Dimension®, Reg. TM: Siemens Healthcare Diagnostic, Glasgow, DE  
 URI-system®, Sure-Vue, Reg. TM: Thermo Fisher Scientific, Waltham, MA  
 Signify®, Contrast, OSOM Reg. TM: Genzyme Corporation, Cambridge, MA  
 KOVA®, Reg. TM: Hycor Biomedical Inc., Garden Grove, CA  
 DiaScreen Series, Reg. TM: Hypoguard, Minneapolis, MN  
 Iris IQ200, Yellow IRIS®, Reg. TM: International Remote Imaging Systems, Inc., Chatsworth, CA  
 Chemstrip®, Chemstrip® 101, Criterion, CUA, Hitachi®, Mical, Mini UA, SUA, URISYS 1100, URISYS Systems, Reg. TM: Roche Diagnostics, Indianapolis, IN  
 QuickVue®, UrinQuick, Reg. TM: Quidel Corporation, San Diego, CA  
 Cen-Slide®, Reg. TM: StatSpin, Norwood, MA  
 PocketChem UA®, Reg. TM: ThermoFisher Scientific, Boulder, CO  
 YD Uriscan®, Reg. TM: YD Diagnostics, Seoul, Korea

\* **Atypical colors may be observed. Results based on intensity.**  
 \* Atypische Farben können beobachtet werden. Resultate basiert worden auf Intensität.  
 \* On peut observer des couleurs atypiques. Résultats basés sur l'intensité.  
 \* Si possono osservare colori atipici. Risultati basati sull'intensità.  
 \* Pueden observarse colores anormales. Los resultados se basan en la intensidad.  
 \* Der kan observeres atypiske farver. Resultater baseret på intensitet.  
 \* Er kan sprake zijn van atypische kleuren. De resultaten zijn gebaseerd op de intensiteit.  
 \* Epätyypillisiä värejä voidaan havaita. Tulokset perustuvat voimakkuuteen.  
 \* Ενδέχεται να παρατηρηθούν άτυπα χρώματα. Τα αποτελέσματα βασίζονται στην ένταση.  
 \* Det kan forekomme atypiske farger. Resultatene er basert på intensitet.  
 \* Može wystąpić niespodziewane zabarwienie. Wyniki na podstawie intensywności.  
 \* Podem ser observadas cores atípicas. Resultados com base na intensidade da cor.  
 \* Наблюдаемый цвет может быть необычным. Результаты зависят от интенсивности.  
 \* Atypiska färger kan observeras. Resultat baserade på intensitet.  
 \* Atipik renkler gözlemlenebilir. Sonuçlar yoğunluğa dayanır.  
 \* 可能观察到非典型颜色。结果依据强度而定。  
 \* 不規則な色が観察される場合があります。結果は明度に基づきます。

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


**If you have any questions concerning this notification, in the USA please contact the Technical Support Department at 800-232-3342 or 510-979-5417. Outside the USA, please contact your local Subsidiary or Distributor.**

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**USA**  
☎ 800-232-3342  
sales.diagnostics.fmt@thermofisher.com

**China**  
☎ +86 800-810-5118  
cdx.cn.info@thermofisher.com

**New Zealand**  
☎ 0800 933 966  
auinfo@thermofisher.com

**Switzerland & Austria**  
☎ +41 26 663 86 70  
cdx.ch.info@thermofisher.com

**Asia Pacific**  
☎ +61 1800 333 110  
cdd.asia.info@thermofisher.com

**France**  
☎ +33 1 40 86 65 20  
cdx.fr.info@thermofisher.com

**Nordic**  
☎ +47 2 325 0433  
info.nordic.cdd@thermofisher.com

**United Kingdom & Ireland**  
☎ +44 1442 868 940  
cdx.uk.info@thermofisher.com

**Australia**  
☎ +61 1800 333 110  
auinfo@thermofisher.com

**Germany**  
☎ +49 0800-40 40 771  
cdx.de.info@thermofisher.com

**South Africa**  
☎ +27117926790  
support-za.idd@thermofisher.com

**Canada**  
☎ 800-282-4075  
info.cddcanada@thermofisher.com

**Japan**  
☎ +81 (0)120-147-075  
JPYOK-CDD.QC@thermofisher.com

**Spain, Portugal & Italy**  
☎ +34 93589 8338  
cdx.es.info@thermofisher.com

**For countries not listed:**  
distributor.cdd@thermofisher.com

**Microgenics Corporation**  
46500 Kato Road  
Fremont, CA 94538 USA  
ISO 13485 Certified Company  
www.thermoscientific.com/diagnostics

**EC REP**

**B.R.A.H.M.S GmbH**  
Neuendorfstrasse 25  
16761 Hennigsdorf, Germany

UAT-INS-VA  
Rev. 14 2016 03

**Thermo**  
SCIENTIFIC