


BAGMATIC LC



Barcode scanner




Stable transport case

Blood donation genius

- Most innovative technology worldwide
- Complete control of all donations
- Programmable error signals
- Data-analysis of donation
- Entirely free adaptable to your needs (volume, speed etc.)
- Free configuration of the software
- Extremely shock resistant

... and many other features!



Technical data

Power Supply	12 V DC
Power Consumption	0,15 A/ 0,35 A
Fuse	T 0.315A ; 250 V
Battery charger	FRIWO FW 1299
Battery charging time	12 hours
Surrounding temperature	10 - 30 ° C
Humidity	maximum 85 %
Dimensions	365 x190x135 mm
Weight	5,1 kg