


Date: 28 February 2006

To whom it may concern:

This is to certify that our Glycohemoglobin (HbA1c) analyzing systems conform to IFCC standard as shown in the chart below.

Traceability chart

IFCC calibrator (1) is assigned using the Primary calibration material and Definitive method by IFCC. The Reference method by ARKRAY is calibrated by IFCC calibrator (1). IFCC value of ARKRAY calibrator (2) is assigned using the Reference method by ARKRAY.

*Reference method by ARKRAY has been reviewed in the IFCC monitoring program and rated as Excellent in the categories of Deviation from IFCC target, Reproducibility, and Linearity.



Tomohide Onishi

Manager

Quality & Safety Management Office

ARKRAY Factory, Inc.

1480 Koji, Konan-cho, Koka, Shiga 520-3306, JAPAN