

Report on Traceability of Standard solution for Urine Sediment Analyzer


Product name: Standard solution for calibration of sediment analyzer

1、Description of traceability

According to ISO 17511 and ICSH, the value of the RBC in standard solution for the urine sediment analyzer are strictly traced to national reference method: single channel, impedance, Z₂ electronic counter.

2、Process of traceability

According to ISO 17511:2003, the following picture shows the tracing chain of the value of RBC in standard solution of the urine sediment analyzer.



Laboratory reference system of Dirui Industrial Co., Ltd.