



**JMS SINGAPORE PTE LTD**  
440 Ang Mo Kio Industrial Park 1, Singapore 569620  
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)  
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)  
Web : www.jmss.com.sg  
Reg : 197901579C

TO MEDISTA SPOL S.R.O.

Page 1 of 1

### Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

CODE NO	DESCRIPTION	LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
811-2577	BLOOD BAG CPDA-1 (S) 250ML IBSP-NP-CLP 16G CE	181003006	2,000	03/10/18	OCT2021
811-8442	BLOOD BAG CPD-SAGM (Q) 450+150+400MLX2 5DPC IBSP-NP-BSB-CLP 16G (DT) CE	180912002	72	12/09/18	SEP2021
811-8442	BLOOD BAG CPD-SAGM (Q) 450+150+400MLX2 5DPC IBSP-NP-BSB-CLP 16G (DT) CE	180925002	3,528	25/09/18	SEP2021
831-8307	BLOOD BAG T&B CPD-SAGM (T) 450+400MLX2 IBSP-NP-BSB-CLP 16G (DT) CE	180911002	3,600	11/09/18	SEP2021

INVOICE NO 1019335  
DATE OF ISSUE 19-OCT-18

TEST ITEMS	JUDGEMENT	REMARKS
1. TOTAL APPEARANCE	PASSED	
2. FUNCTION TEST	PASSED	
3. PYROGENITY-USP	PASSED	
4. STERILITY-USP	PASSED	
5. SAFETY/TOXICITY-USP	PASSED	
6. CHEMICAL CONTENTS	PASSED	

TEST ITEMS 3-6 ARE NOT APPLICABLE TO BULK AND NON-STERILE PRODUCTS  
TEST ITEM 6 IS APPLICABLE TO BLOOD BAG ONLY

P.O. NO : 2180299,2180328

CY CHIA  
QA MANAGER



JMS SINGAPORE PTE LTD  
440 Ang Mo Kio Industrial Park 1, Singapore 569620  
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)  
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)  
Web : www.jmss.com.sg  
Reg : 197901579C

PRODUCT CODE : 811-8442  
LOT NUMBER : 180925002  
DESCRIPTION : BLOOD BAG CPD-SAGM (Q) 450+150+400MLX2 5DPC IBSP-NP-BSB-CLP 16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.4	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	100.6	PASSED
	* TOTAL CITRATE (%)	95 - 105	102.2	PASSED
	DEXTROSE (%)	95 - 105	100.7	PASSED
	TOTAL SODIUM (%)	95 - 105	99.7	PASSED

\* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.5	5.3	PASSED
	SODIUM (%)	95 - 105	100.0	PASSED
	CHLORIDE (%)	95 - 105	101.3	PASSED
	DEXTROSE (%)	95 - 105	100.5	PASSED
	MANNITOL (%)	95 - 105	100.2	PASSED
	ADENINE (%)	90 - 110	100.0	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 19-OCT-18

CY CHIA  
QA MANAGER