

JMSS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

Page 1 of 1

TO MEDISTA SPOL S.R.O.

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS


INVOICE NO 1019971
DATE OF ISSUE 13-MAY-19

CODE NO	DESCRIPTION	LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
811-4032	BLOOD BAG CPDA-1 (S) 450ML IBSP-NP-CLP 16G CE	190409003	240	09/04/19	APR2022
814-0132	TRANSFER BAG 150ML CE	190301002	200	01/03/19	MAR2022

TEST ITEMS	JUDGEMENT	REMARKS
1. TOTAL APPEARANCE	PASSED	
2. FUNCTION TEST	PASSED	
3. PYROGENITY-USP	PASSED	
4. STERILITY-USP	PASSED	
5. SAFETY/TOXICITY-USP	PASSED	
6. CHEMICAL CONTENTS	PASSED	

P.O. NO : 2190018,2190065

TEST ITEMS 3-6 ARE NOT APPLICABLE TO BULK AND NON-STERILE PRODUCTS
TEST ITEM 6 IS APPLICABLE TO BLOOD BAG ONLY


CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
 440 Ang Mo Kio Industrial Park 1, Singapore 569620
 Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
 Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
 Web : www.jmss.com.sg
 Reg : 197901579C

PRODUCT CODE : 811-4032

LOT NUMBER : 190409003

DESCRIPTION : BLOOD BAG CPDA-1 (S) 450ML IBSP-NP-CLP 16G CE

TEST ITEMS		SPECIFICATION LIMITS	RESULTS
CHEMICAL CONTENTS	pH	5.0 - 6.0	PASSED 5.7
	CITRIC ACID (%)	95 - 105	PASSED 99.4
	SODIUM BIPHOSPHATE (%)	95 - 105	PASSED 100.4
	* TOTAL CITRATE (%)	95 - 105	PASSED 101.1
	DEXTRROSE (%)	95 - 105	PASSED 99.6
	ADENINE (%)	90 - 110	PASSED 99.5

* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 13-MAY-19

CY CHIA

QA MANAGER