



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

TO BIONIC MEDIZINTECHNIK G.M.B.H.

Page 1 of 4

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

INVOICE NO 1020488, 490
DATE OF ISSUE 08-OCT-19

CODE NO	DESCRIPTION	LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
420-1114	DONOR NEEDLE 16G NON-STERILE	190620535	12,000	20/06/19	-
420-1114	DONOR NEEDLE 16G NON-STERILE	190621531	24,000	21/06/19	-
420-1114	DONOR NEEDLE 16G NON-STERILE	190621532	24,000	21/06/19	-
420-1114	DONOR NEEDLE 16G NON-STERILE	190624531	24,000	24/06/19	-
420-1114	DONOR NEEDLE 16G NON-STERILE	190624532	24,000	24/06/19	-
420-1114	DONOR NEEDLE 16G NON-STERILE	190625531	24,000	25/06/19	-
420-1114	DONOR NEEDLE 16G NON-STERILE	190625532	24,000	25/06/19	-
420-1114	DONOR NEEDLE 16G NON-STERILE	190626531	12,000	26/06/19	-
811-1000	CORD BLOOD COLLECTION SET CPDA-1 (S) 157ML W/INJ PORT NP-CLP 16G CE	190913001	3,150	13/09/19	SEP2022
814-0302	TRANSFER BAG 300ML CE	190903001	4,000	03/09/19	SEP2022
814-0302	TRANSFER BAG 300ML CE	190911001	2,000	11/09/19	SEP2022
822-1616	SYSLOC AVF 16GX1" BE 15CM W/CLAMP CE	190822211	36,000	22/08/19	AUG2024



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

TO BIONIC MEDIZINTECHNIK G.M.B.H.

Certificate of Compliance

Page 2 of 4

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

INVOICE NO 1020488, 490						DATE OF ISSUE 08-OCT-19			
DESCRIPTION						LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
822-1616	SYSLOC AVF 16GX1" BE 15CM W/CLAMP CE					190822212	36,000	22/08/19	AUG2024
822-1616	SYSLOC AVF 16GX1" BE 15CM W/CLAMP CE					190823211	36,000	23/08/19	AUG2024
822-1616	SYSLOC AVF 16GX1" BE 15CM W/CLAMP CE					190828212	36,000	28/08/19	AUG2024
822-1616	SYSLOC AVF 16GX1" BE 15CM W/CLAMP CE					190830211	36,000	30/08/19	AUG2024
822-1616	SYSLOC AVF 16GX1" BE 15CM W/CLAMP CE					190830212	36,000	30/08/19	AUG2024
822-1616	SYSLOC AVF 16GX1" BE 15CM W/CLAMP CE					190912211	36,000	12/09/19	SEP2024
822-1616	SYSLOC AVF 16GX1" BE 15CM W/CLAMP CE					190913211	36,000	13/09/19	SEP2024
862-1506	SYSLOC AVF WITH TERMINAL LOCK 15GX4/5" BE 15CM W/CLAMP CE					190913212	36,000	13/09/19	SEP2024
862-1507	SYSLOC AVF WITH TERMINAL LOCK 15GX4/5" 15CM W/CLAMP CE					190911212	6,000	11/09/19	SEP2024
862-1510	SYSLOC AVF WITH TERMINAL LOCK 15GX4/5" BE 30CM W/CLAMP CE					190911214	7,000	11/09/19	SEP2024
862-1513	SYSLOC AVF WITH TERMINAL LOCK 15GX4/5" 30CM W/CLAMP CE					190904231	4,000	04/09/19	SEP2024
						190904232	8,000	04/09/19	SEP2024



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64595564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

TO BIONIC MEDIZINTECHNIK G.M.B.H.


Page 3 of 4

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

INVOICE NO 1020488, 490
DATE OF ISSUE 08-OCT-19

CODE NO	DESCRIPTION	LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
862-1515	SYSLOC AVF WITH TERMINAL LOCK 15GX1" 30CM W/CLAMP CE	190905231	1,000	05/09/19	SEP2024
862-1518	SYSLOC AVF WITH TERMINAL LOCK 15GX4/5" 50CM W/CLAMP CE	190913332	1,200	13/09/19	SEP2024
862-1609	SYSLOC AVF WITH TERMINAL LOCK 16GX4/5" BE 30CM W/CLAMP CE	190904233	9,000	04/09/19	SEP2024
862-1610	SYSLOC AVF WITH TERMINAL LOCK 16GX4/5" 30CM W/CLAMP CE	190904234	1,000	04/09/19	SEP2024
862-1613	SYSLOC AVF WITH TERMINAL LOCK 16GX1" BE 15CM W/CLAMP CE	190910211	2,000	10/09/19	SEP2024
862-1614	SYSLOC AVF WITH TERMINAL LOCK 16GX1" BE 30CM W/CLAMP CE	190905232	8,000	05/09/19	SEP2024
862-1616	SYSLOC AVF WITH TERMINAL LOCK 16GX1" 30CM W/CLAMP CE	190905233	4,000	05/09/19	SEP2024
862-1617	SYSLOC AVF WITH TERMINAL LOCK 16GX3/5" BE 30CM W/CLAMP CE	190905234	1,000	05/09/19	SEP2024
862-1700	SYSLOC AVF WITH TERMINAL LOCK 17GX4/5" BE 15CM W/CLAMP CE	190911211	2,000	11/09/19	SEP2024
862-1707	SYSLOC AVF WITH TERMINAL LOCK 17GX4/5" BE 30CM W/CLAMP CE	190916231	3,000	16/09/19	SEP2024
862-1708	SYSLOC AVF WITH TERMINAL LOCK 17GX4/5" 30CM W/CLAMP CE	190916232	2,000	16/09/19	SEP2024
862-1711	SYSLOC AVF WITH TERMINAL LOCK 17GX1" BE 15CM W/CLAMP CE	190911213	1,000	11/09/19	SEP2024



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

Page 4 of 4

TO BIONIC MEDIZINTECHNIK G.M.B.H.

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

TEST ITEMS	JUDGEMENT	REMARKS
1. TOTAL APPEARANCE	PASSED	
2. FUNCTION TEST	PASSED	
3. PYROGENITY-USP	PASSED	
4. STERILITY-USP	PASSED	
5. SAFETY/TOXICITY-USP	PASSED	
6. CHEMICAL CONTENTS	PASSED	

TEST ITEMS 3-6 ARE NOT APPLICABLE TO BULK AND NON-STERILE PRODUCTS
TEST ITEM 6 IS APPLICABLE TO BLOOD BAG ONLY

INVOICE NO 1020488, 490
DATE OF ISSUE 08-OCT-19

P.O. NO : 2190276,2190263,2190231,2190292

CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

PRODUCT CODE : 811-1000

LOT NUMBER : 190913001

DESCRIPTION : CORD BLOOD COLLECTION SET CPDA-1 (S) 157ML W/INJ PORT NP-CLP 16G CE

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.4	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	100.7	PASSED
	* TOTAL CITRATE (%)	95 - 105	101.4	PASSED
	DEXTROSE (%)	95 - 105	99.4	PASSED
	ADENINE (%)	90 - 110	99.0	PASSED

* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 08-OCT-19

CY CHIA

QA MANAGER