



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

Page 1 of 2

TO MEDISTA SPOL S.R.O.

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

INVOICE NO 1021183
DATE OF ISSUE 31-MAR-20

CODE NO	DESCRIPTION	LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
201-0111	BB PLATE CLAMP NO.2 JMS-S NON-STERILE	191117161	3,600	17/11/19	-
811-8442	BLOOD BAG CPD-SAGM (Q) 450+150+400MLX2 5DPC IBSP-NP-BSB-CLP 16G (DT) CE	200213002	3,600	13/02/20	FEB2023
814-0503	TRANSFER BAG 500ML CE	200109001	1,000	09/01/20	JAN2023
814-0504	TRANSFER BAG 500ML W/Y CONNECTOR CE	200110001	1,040	10/01/20	JAN2023
831-8307	BLOOD BAG T&B CPD-SAGM (T) 450+400MLX2 IBSP-NP-BSB-CLP 16G (DT) CE	200117002	2,640	17/01/20	JAN2023
833-8441	BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF 16G (DT) CE	200212001	396	12/02/20	FEB2023
833-8441	BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF 16G (DT) CE	200226002	468	26/02/20	FEB2023
833-8441	BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF 16G (DT) CE	200304002	324	04/03/20	MAR2023



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

Page 2 of 2

TO MEDISTA SPOL S.R.O.

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

TEST ITEMS	JUDGEMENT	REMARKS
1. TOTAL APPEARANCE	PASSED	
2. FUNCTION TEST	PASSED	
3. PYROGENITY-USP	PASSED	
4. STERILITY-USP	PASSED	
5. SAFETY/TOXICITY-USP	PASSED	
6. CHEMICAL CONTENTS	PASSED	

TEST ITEMS 3-6 ARE NOT APPLICABLE TO BULK AND NON-STERILE PRODUCTS
TEST ITEM 6 IS APPLICABLE TO BLOOD BAG ONLY

INVOICE NO 1021183 **DATE OF ISSUE** 31-MAR-20

P.O. NO : 2190385,2190349,2190424

CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

PRODUCT CODE : 811-8442
LOT NUMBER : 200213002
DESCRIPTION : BLOOD BAG CPD-SAGM (Q) 450+150+400MLX2 5DPC IBSP-NP-BSB-CLP 16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.4	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	99.2	PASSED
	* TOTAL CITRATE (%)	95 - 105	100.4	PASSED
	DEXTROSE (%)	95 - 105	98.5	PASSED
	TOTAL SODIUM (%)	95 - 105	99.5	PASSED

* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	4.5 - 6.5	5.2	PASSED
	SODIUM (%)	95 - 105	99.5	PASSED
	CHLORIDE (%)	95 - 105	100.6	PASSED
	DEXTROSE (%)	95 - 105	100.0	PASSED
	MANNITOL (%)	95 - 105	100.1	PASSED
	ADENINE (%)	90 - 110	101.1	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 31-MAR-20

CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

PRODUCT CODE : 831-8307

LOT NUMBER : 200117002

DESCRIPTION : BLOOD BAG T&B CPD-SAGM (T) 450+400MLX2 IBSP-NP-BSB-CLP 16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.3	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	100.1	PASSED
	* TOTAL CITRATE (%)	95 - 105	100.8	PASSED
	DEXTROSE (%)	95 - 105	99.0	PASSED
	TOTAL SODIUM (%)	95 - 105	99.2	PASSED

* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	4.5 - 6.5	5.2	PASSED
	SODIUM (%)	95 - 105	99.6	PASSED
	CHLORIDE (%)	95 - 105	100.9	PASSED
	DEXTROSE (%)	95 - 105	99.3	PASSED
	MANNITOL (%)	95 - 105	102.1	PASSED
	ADENINE (%)	90 - 110	99.6	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 31-MAR-20

CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

PRODUCT CODE : 833-8441

LOT NUMBER : 200212001

DESCRIPTION : BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF
16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.4	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	98.8	PASSED
	* TOTAL CITRATE (%)	95 - 105	100.6	PASSED
	DEXTROSE (%)	95 - 105	98.9	PASSED
	TOTAL SODIUM (%)	95 - 105	99.3	PASSED

* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	4.5 - 6.5	5.2	PASSED
	SODIUM (%)	95 - 105	99.9	PASSED
	CHLORIDE (%)	95 - 105	100.2	PASSED
	DEXTROSE (%)	95 - 105	99.6	PASSED
	MANNITOL (%)	95 - 105	100.4	PASSED
	ADENINE (%)	90 - 110	100.5	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 31-MAR-20

CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

PRODUCT CODE : 833-8441
LOT NUMBER : 200226002
DESCRIPTION : BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF
16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.4	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	99.7	PASSED
	* TOTAL CITRATE (%)	95 - 105	100.6	PASSED
	DEXTROSE (%)	95 - 105	98.6	PASSED
	TOTAL SODIUM (%)	95 - 105	100.7	PASSED


* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	4.5 - 6.5	5.2	PASSED
	SODIUM (%)	95 - 105	100.8	PASSED
	CHLORIDE (%)	95 - 105	100.6	PASSED
	DEXTROSE (%)	95 - 105	100.5	PASSED
	MANNITOL (%)	95 - 105	100.3	PASSED
	ADENINE (%)	90 - 110	100.9	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 31-MAR-20



CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

PRODUCT CODE : 833-8441

LOT NUMBER : 200304002

DESCRIPTION : BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF
16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.3	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	100.0	PASSED
	* TOTAL CITRATE (%)	95 - 105	100.6	PASSED
	DEXTROSE (%)	95 - 105	98.8	PASSED
	TOTAL SODIUM (%)	95 - 105	100.1	PASSED

* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	4.5 - 6.5	5.2	PASSED
	SODIUM (%)	95 - 105	100.2	PASSED
	CHLORIDE (%)	95 - 105	100.5	PASSED
	DEXTROSE (%)	95 - 105	100.4	PASSED
	MANNITOL (%)	95 - 105	100.8	PASSED
	ADENINE (%)	90 - 110	100.8	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 31-MAR-20

CY CHIA
QA MANAGER