



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

Page 1 of 2

TO BIONIC MEDIZINTECHNIK G.M.B.H.

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

INVOICE NO DATE OF ISSUE
1022399 17-FEB-21

CODE NO	DESCRIPTION	LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
811-1000	CORD BLOOD COLLECTION SET CPDA-1 (S) 157ML W/INJ PORT NP-CLP 16G CE	201109004	3,150	09/11/20	NOV2023
814-0107	TRANSFER BAG 1000ML CE	210104001	2,000	04/01/21	JAN2024
814-0121	TRANSFER BAG 150ML W/4 BAGS CE	201229001	3,000	29/12/20	DEC2023
814-0121	TRANSFER BAG 150ML W/4 BAGS CE	201231001	3,000	31/12/20	DEC2023
814-0132	TRANSFER BAG 150ML CE	210105001	2,000	05/01/21	JAN2024
814-0302	TRANSFER BAG 300ML CE	210111001	1,000	11/01/21	JAN2024
814-0302	TRANSFER BAG 300ML CE	210122001	9,000	22/01/21	JAN2024
814-0611	TRANSFER BAG 600ML CE	210108001	2,000	08/01/21	JAN2024
814-0611	TRANSFER BAG 600ML CE	210127001	3,000	27/01/21	JAN2024
831-8411	BLOOD BAG T&B CPD-SAGM (Q) 450+400MLX3 5DPC IBSP-NP-BSB-CLP 16G (DT) CE	210113002	1,000	13/01/21	JAN2024
833-8441	BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF 16G (DT) CE	201229001	1,260	29/12/20	DEC2023



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

Page 2 of 2

TO BIONIC MEDIZINTECHNIK G.M.B.H.

Certificate of Compliance

WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS

INVOICE NO DATE OF ISSUE
1022399 17-FEB-21

P.O. NO : 2200400,2200322,2200397,2200
364,2200362,2200320,2200348

TEST ITEMS	JUDGEMENT	REMARKS
1.TOTAL APPEARANCE	PASSED	
2.FUNCTION TEST	PASSED	
3.PYROGENITY-USP	PASSED	
4.STERILITY-USP	PASSED	
5. SAFETY/TOXICITY-USP	PASSED	
6. CHEMICAL CONTENTS	PASSED	

TEST ITEMS 3-6 ARE NOT APPLICABLE TO BULK AND NON-STERILE PRODUCTS
TEST ITEM 6 IS APPLICABLE TO BLOOD BAG ONLY

CY CHIA
QA MANAGER



JMS SINGAPORE PTE LTD
440 Ang Mo Kio Industrial Park 1, Singapore 569620
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)
Web : www.jmss.com.sg
Reg : 197901579C

PRODUCT CODE : 833-8441

LOT NUMBER : 201229001

DESCRIPTION : BLOOD BAG T&T CPD-SAGM (Q) 450+450+400+450ML IBSP-BSB-NP-CLP-WBF
16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.7	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	99.8	PASSED
	* TOTAL CITRATE (%)	95 - 105	101.4	PASSED
	DEXTROSE (%)	95 - 105	99.5	PASSED
	TOTAL SODIUM (%)	95 - 105	100.3	PASSED

* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	4.5 - 6.5	5.3	PASSED
	SODIUM (%)	95 - 105	100.7	PASSED
	CHLORIDE (%)	95 - 105	100.7	PASSED
	DEXTROSE (%)	95 - 105	101.5	PASSED
	MANNITOL (%)	95 - 105	99.5	PASSED
	ADENINE (%)	90 - 110	101.4	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 17-FEB-21

CY CHIA
QA MANAGER