



**JMS SINGAPORE PTE LTD**  
440 Ang Mo Kio Industrial Park 1, Singapore 569620  
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)  
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)  
Web : www.jmss.com.sg  
Reg : 197901579C

Page 1 of 1

**TO MEDISTA SPOL S.R.O.**

**Certificate of Compliance**

**WE CERTIFY ALL LOTS MEET JMSS SPECIFICATIONS**

**INVOICE NO** 1023355  
**DATE OF ISSUE** 13-DEC-21

CODE NO	DESCRIPTION	LOT NO	QUANTITY	PROD. DATE	EXPIRY DATE
811-8442	BLOOD BAG CPD-SAGM (Q) 450+150+400MLX2 5DPC IBSP-NP-BSB-CLP 16G (DT) CE	211021003	2,052	21/10/21	OCT2024
814-0132	TRANSFER BAG 150ML CE	210923001	1,000	23/09/21	SEP2024
814-0503	TRANSFER BAG 500ML CE	211001001	1,000	01/10/21	OCT2024

TEST ITEMS	JUDGEMENT	REMARKS
1.TOTAL APPEARANCE	PASSED	
2.FUNCTION TEST	PASSED	
3.PYROGENITY-USP	PASSED	
4.STERILITY-USP	PASSED	
5. SAFETY/TOXICITY-USP	PASSED	
6. CHEMICAL CONTENTS	PASSED	

TEST ITEMS 3-6 ARE NOT APPLICABLE TO BULK AND NON-STERILE PRODUCTS  
TEST ITEM 6 IS APPLICABLE TO BLOOD BAG ONLY

**P.O. NO** : 2210231,2210268

CY CHIA  
QA MANAGER



JMS SINGAPORE PTE LTD  
440 Ang Mo Kio Industrial Park 1, Singapore 569620  
Tel : (65) 64571176 (Main), (65) 64571144 (Sales)  
Fax : (65) 64548500 (Main), (65) 64599564 (Sales)  
Web : www.jmss.com.sg  
Reg : 197901579C

PRODUCT CODE : 811-8442  
LOT NUMBER : 211021003  
DESCRIPTION : BLOOD BAG CPD-SAGM (Q) 450+150+400MLX2 5DPC IBSP-NP-BSB-CLP 16G (DT) CE

CPD

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	5.0 - 6.0	5.7	PASSED
	CITRIC ACID (%)	95 - 105	99.7	PASSED
	SODIUM BIPHOSPHATE (%)	95 - 105	99.6	PASSED
	* TOTAL CITRATE (%)	95 - 105	100.9	PASSED
	DEXTROSE (%)	95 - 105	100.8	PASSED
	TOTAL SODIUM (%)	95 - 105	99.3	PASSED

\* TOTAL CITRATE = CITRIC ACID + SODIUM CITRATE

SAGM

TEST ITEMS		SPECIFICATION LIMITS	RESULTS	
CHEMICAL CONTENTS	pH (%)	4.5 - 6.5	5.2	PASSED
	SODIUM (%)	95 - 105	100.1	PASSED
	CHLORIDE (%)	95 - 105	100.7	PASSED
	DEXTROSE (%)	95 - 105	99.2	PASSED
	MANNITOL (%)	95 - 105	100.0	PASSED
	ADENINE (%)	90 - 110	101.1	PASSED

THE ABOVE RESULTS ARE WITHIN SPECIFICATIONS.

DATE OF ISSUE : 13-DEC-21

CY CHIA  
QA MANAGER