

TERUMOBCT

Supplier:

Terumo BCT, Inc.
10811 West Collins Avenue
Lakewood, CO 80215, USA
Toll Free (US & Canada): 877-339-4228
Phone :303-231-4357
FAX :303-542-5215

Manufacturer:

TERUMO CORPORATION Fujinomiya Factory
818, Misonodaira, Fujinomiya-shi, Shizuoka 418-0004 JAPAN
Phone :+81-544-27-7142
FAX :+81-544-23-4548
Compliance and Quality Administration Department:
Yoshiki Takagi

LOT CERTIFICATE

Product Name:	TERUFLEX TRANSFER BAG		
Product Code/Cat No:	BB*T100BM		
Lot number:	190410F3		
Lot size (amount):	3,200 sets		
Sterilization Date:	11/04/2019	Sterilization Method:	Ethylene Oxide
Expiry:	03/2022		

Material validation	Quality of product is validated according to ISO3826-1.
----------------------------	---

QC specifications and inspection result	(1)Chemical test for PVC material based on internal standards	<input checked="" type="checkbox"/>	Passed
	(2)Chemical test for other plastic material based on internal standards	<input checked="" type="checkbox"/>	Passed
	(3)Function test for PVC container based on internal standards	<input checked="" type="checkbox"/>	Passed
	(4)Biological test (Bacterial Endotoxin) according to ISO 3826-1	<input checked="" type="checkbox"/>	Passed
	(5)Microbiological test (Biological Indicator test) according to EN ISO 11135	<input checked="" type="checkbox"/>	Passed

Comments results: Lot No. 190410F3 meets all requirements of QC product specification
Product Code/Cat No: BB*T100BM

Release statement : Lot Released

Date of release determination : 09 / 05 / 2019
(Day) (Month) (Year)



Date of signature : 28 / 05 / 2019
(Day) (Month) (Year)

Yoshiki Takagi,
Compliance and Quality Administration Department
Terumo Corporation Fujinomiya Factory

Terumo BCT, Inc. USA
10811 West Collins Ave., Lakewood, CO 80215-4440
USA Phone: +1.877.339.4228
Phone: +1.303.231.4357, Fax: +1.303.542.5215
TERUMOBCT.COM

Terumo BCT Europe N.V.
Europe, Middle East and Africa
Ikaroslaan 41
1930 Zaventem, Belgium
Phone: +32.2.715.05.90, Fax: +32.2.721.07.70

Terumo BCT (Asia Pacific) Ltd.
Room 3903-3903A, 39/F
New York Life Tower, Windsor House
311 Gloucester Road, Causeway Bay, Hong Kong
Phone: +852.2283.0700, Fax: +852.2576.1311