

Biosealer CR 6




Biosealer CR 6

Biosealer CR 6 is a battery operated sealing unit for PVC-tubes, especially bloodbags tubes or sets for plasmaferes. With new and fully loaded batteries it manages more than 1000 seals on PVC-tubes. The high capacity makes it suitable to use CR 6 both in donation rooms, laboratories and in mobile use. Biosealer CR 6 is comprised of 3 units:

- Sealing handle with coax cable 1,9 m.
- Power unit with batteries, Rf generator and a control unit monitored by a micro processor.
- Battery charger.

Optional there is a carry bag for mobile use.

Technical data:

Batteries:	2 x 12V 2Ah
Internal fuse:	Manufacture: Sonnenschein A500
Effect consumption:	Overload protection type PTC, self recovery
Sealing time:	120 W maximum effect
Frequency:	Adjustable 0,5 - 5 sec
Output effect:	40,68 MHz crystal controlled
Dimensions:	70 W maximum effect
Weight:	266 x 119 x 88 mm (L x W x H) incl. rf-outlet
Classification:	3,1 kg incl. batteries
Safety approval:	Protection against electrical chock:
EMC-approval:	class II type B. 
	EN 60601-1 (IEC601-1)
	EN 60601-1-2 (IEC601-1-2),
	emission EN55011 grupp2 class B

LJUNGBERG & KÖGEL AB

Box 1032, 251 10 Helsingborg, Sweden